

CFPM100

IN STOCK PRODUCT
 Visit www.oakworksmed.com for pricing.



- Best Adjustable Height Range **22"- 40" Industry Lowest**
- 56" Metal Free Imaging Space
- 500 lb. Patient Weight Capacity
- **4 Year Warranty**

PRODUCT SPECS

PATIENT WEIGHT CAPACITY	500 lbs. (227kg.)
INTEGRATED HEAD REST CARBON FIBER TOP	22" W x 84" L (56cm x 213cm)
MOTION height range 22" - 40" (56cm x 102cm)	
CONTROLS	1 foot
TABLE WEIGHT	350 lbs. (170kg.)
WARRANTY	4 year
REGULATORY	FDA, ETL, CE

The OAKWORKS® CFPM100 C-Arm Imaging Table uses our robust single tower lift system to create this 500 lbs patient weight capacity table. This cost effective 22" x 84" carbon fiber table is ideal for large patient care in pain management practices, surgical centers, and other multi-disciplinary medical facilities that do not require 2, 3 or 4 table movements. The 22" - 40" adjustable electric height range is the lowest in the mobile imaging table market and promotes safe and secure wheelchair transfer. Accessories included: 1" radiolucent table top pad (2" pad is optional), 1 crescent face pad, 2 supine cranial pads and 1 padded safety strap. See page 1 for pad and strap descriptions.



ADJUSTABLE HEIGHT RANGE
 Table height range of 22"- 40"
INDUSTRY'S LOWEST!



HAND RAIL
 Move your table with ease

RECTANGULAR TOP OPTION

The rectangular top is wider and shorter (24"x80", 61cm-203 cm) than the integrated headrest top providing greater flexibility.

The three cranial positioning pads are not included with this option.



Made in the USA 
 with US & imported parts

OAKWORKS®
MEDICAL

C-ARM IMAGING TABLES