

CFPM400

4 Movement Table

IN STOCK PRODUCT
 Visit www.oakworksmed.com for pricing.

- **Best Adjustable Height Range 26"- 44" Industry's Lowest**
- 56" Metal Free Imaging Space
- Cost Effective 4 Movement Table
- 500 lbs. Patient Weight Capacity
- **4 Year Warranty**



view online video
www.oakworksmed.com

PRODUCT SPECS

PATIENT WEIGHT CAPACITY	500 lb. (227kg.)
CARBON FIBER TOP	22" x 84" (56cm x 213cm)
MOTION 1	height range 26" - 44" (66cm x 112cm)
MOTION 2	10" overall travel top (12cm)
MOTION 3	Trendelenburg tilt +/- 12°/15°
MOTION 4	lateral tilt +/- 15°
CONTROLS	1 foot, 1 hand
TABLE WEIGHT	375 lbs. (170 kg.)
WARRANTY	4 year
REGULATORY	FDA, ETL, CE

see table features on next page



INCLUDED ACCESSORIES:

Patient Safety Strap
 Adjustable Velcro® strap to help immobilize the patient during anesthesia and surgery.

1" Radiolucent Table top pad (2" pad is optional)

3 Cranial Positioning Pads

The three pads provide versatile support in supine and prone positions.

- 2" (5 cm) thick x 12" (31 cm) (square)
- 4" (10 cm) thick x 12" (31 cm) (square)
- Crescent face rest pad



Made in the USA 
 with US & imported parts

OAKWORKS®
MEDICAL

C-ARM IMAGING TABLES

CFPM400

4 Movement Table Features

The OAKWORKS® CFPM400 C-Arm Imaging Table is the best valued 4 movement imaging table in the mobile imaging table market. Incorporating OAKWORKS® brand reputation for quality & reliability, the CFPM400 is designed to accommodate mobile c-arm and ceiling suspended imaging systems. It's dynamic height range of 26" to 44" accommodates low height patient wheelchair transfers to maximum height c-arm requirements for full lateral imaging. Mounted on a low profile base, this mobile table can be secured via 4 locking casters and provides 56" of metal free imaging area. The CFPM400 is FDA listed, ETL approved and CE marked.



HEIGHT RANGE 26"- 44"



TRENDELENBURG TILT +12°/-15°



TRAVEL TOP 10"



LATERAL TILT +/-15°

RECTANGULAR TOP OPTION

The rectangular top is wider and shorter (24"x80", 61cm-203 cm) than the integrated headrest top providing greater flexibility to customize the configuration of our imaging tables to better meet the patient positioning needs of any medical imaging facility.

The three cranial positioning pads are not included with this option.

