

# DTPM300

## 3 Movement Table

IN STOCK PRODUCT  
 Visit [www.oakworksmed.com](http://www.oakworksmed.com) for pricing.



### PRODUCT SPECS

PATIENT WEIGHT CAPACITY	450lbs. (205kg.)
TOP	24" (61 cm) W, 78" (198 cm) L
MOTION 1	height range 24" - 40" (61 cm x 102cm)
MOTION 2	23" (56 cm) overall manual travel top
MOTION 3	Trendelenburg tilt +/-14°
CONTROLS	1 foot, 1 hand
TABLE WEIGHT	320lbs. (813kg.)
WARRANTY	4 year
REGULATORY	FDA, ETL, CE

- Most Cost-effective 3 Movement Imaging Table
- 2 Tower Lift Design
- 40" Open Space Between Towers
- 4 Year Warranty

OAKWORKS® DTPM300 is ideal for pain management treatment centers, surgical centers, and other multi-disciplinary medical facilities where a cost effective solution without compromise is needed. Accessories included: 1" radiolucent table top pad (2" pad is optional)

### MOVEMENTS



Made in the USA   
 with US & imported parts

OAKWORKS®  
**MEDICAL**

C-ARM IMAGING TABLES